

100mM dNTP Set

Shipping: Dry Ice/Blue Ice	Catalog numbers	
Batch No.: See vial	BIO-39025: 4 x 25 µmol	4 x 250 µl
Concentration: 100 mM each	BIO-39026: 4 x 100 µmol	4 x 4 x 250 µl
	BIO-39049: 4 x 100 µmol	4 x 1 ml



A Meridian Life Science® Company

Store at -20°C

Storage and stability:

dNTP Set is shipped on dry/blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided. For long-term storage, aliquoting is recommended.

Expiry:

When stored under the recommended conditions and handled correctly, quality is retained until the expiry date on the outer box label.

Safety precautions:

Please refer to the material safety data sheet for information.

Presentation:

A set of 4 separate 100 mM lithium salt solutions (dATP, dGTP, dCTP, and dTTP, (pH 7.5). Each solution contains 25 µmol (250 µl) of the corresponding dNTP.

Quality control specifications:

dNTP Mix is extensively tested for contamination prior to release. dNTPs are >99% pure by HPLC and are free of DNase, RNase, Protease, phosphatase and nicking activity. Lithium salts, >99% deoxynucleoside triphosphates (HPLC, area %), <1% deoxynucleoside mono-phosphates and deoxynucleoside diphosphates.

Description

A set of ready-to-use molecular grade dNTP solutions consisting of 4 separate 100mM solutions of dATP, dGTP, dCTP, and dTTP, at pH 7.5 and supplied as lithium salts in purified water. For use in DNA polymerization reactions, DNA labeling and sequencing processes. Dependable PCR grade.

All Bioline dNTPs are supplied as lithium salts in purified water at pH 7.5. Lithium salts have greater resistance to repeated freezing and thawing cycles than sodium salts, and lithium salt dNTP preparations remain sterile over the entire shelf-life due to the bacteriostatic activity of lithium towards various microorganisms.

This reagent has been manufactured under 13485 Quality Management System and is suitable for further manufacturing use as an IVD component.

Features

- Ultra-pure: >99% trisphosphate by HPLC
- Extended shelf-life of 24 months at -20 °C
- Free from PCR inhibitors
- DNase, RNase and Nickase free
- Manufactured by Bioline in a purpose-built facility

Applications

- Standard and long range PCR assays
- cDNA synthesis
- qPCR
- Microarrays
- DNA sequencing
- Labeling

dNTP Master Mix Preparation

Prepare a Master Mix for DNA synthesis as follows: Mix equal volumes of all separate dNTP solutions in a new micro-centrifuge tube. The final solution has a concentration of 25mM of each dNTP, which corresponds to a 100x working concentration.

Product citations

1. Rip, D., *et al. Food Analyt. Meth.* **2**, 190-196 (2009).
2. Thomson, S., *et al. Meth. Mol. Biol.* **512**, 233-248 (2009).
3. Hayden, M.J., *et al. BMC Genom.* **9**,80-85 (2008).
4. Mayor, N.P., *et al. J. Immunol. Meth.* **327(1/2)**, 82-87 (2007).

Associated products

Product Name	Pack Size	Catalog No.
HyperLadder™ 1 Kb	200 Lanes	BIO-33025
VELOCITY DNA Polymerase	250 Units	BIO-21098
Agarose	500g	BIO-41025

	dATP	dCTP	dGTP	dTTP
Product	dATP Lithium 100mM Solution	dCTP Lithium 100mM Solution	dGTP Lithium 100mM Solution	dTTP Lithium 100mM Solution
Nomenclature	2'-deoxyadenosine-5'-triphosphate	2'-deoxycytidine-5'-triphosphate	2'-deoxyguanosine-5'-triphosphate	2'-deoxythymidine-5'-triphosphate
Formula	C ₁₀ H ₁₂ N ₅ O ₁₂ P ₃ Li ₄	C ₉ H ₁₂ N ₃ O ₁₃ P ₃ Li ₄	C ₁₀ H ₁₂ N ₅ O ₁₃ P ₃ Li ₄	C ₁₀ H ₁₃ N ₂ O ₁₄ P ₃ Li ₄
Molecular Weight	514.9 g/mol	490.9 g/mol	530.9 g/mol	505.9 g/mol
λ _{max} pH 7.0	259nm	272nm	252nm	267nm
ε at λ _{max} @ pH7.0	15.4 E x mmol ⁻¹ x cm ⁻¹	9.1 E x mmol ⁻¹ x cm ⁻¹	13.7 E x mmol ⁻¹ x cm ⁻¹	9.6 E x mmol ⁻¹ x cm ⁻¹
A ₂₅₀ /A ₂₆₀	0.78 ± 0.03	0.82 ± 0.03	1.16 ± 0.05	0.65 ± 0.03
A ₂₈₀ /A ₂₆₀	0.15 ± 0.02	0.98 ± 0.03	0.66 ± 0.03	0.73 ± 0.02
Concentration	100 mM ± 2%	100 mM ± 2%	100 mM ± 2%	100 mM ± 2%
Appearance	Clear Colorless Solution	Clear Colorless Solution	Clear Colorless Solution	Clear Colorless Solution
pH of Solution	7.5	7.5	7.5	7.5
dNTP (HPLC Area)	≥99%	≥99%	≥99%	≥99%
dNDP (HPLC Area)	<1%	<1%	<1%	<1%
DNases, RNases, Nicking Activity	Negative	Negative	Negative	Negative
Storage	at -20 °C	at -20 °C	at -20 °C	at -20 °C
Stability	≤24 months	≤24 months	≤24 months	≤24 months

Bioline Reagents Ltd
UNITED KINGDOM

Tel: +44 (0)20 8830 5300

Bioline USA Inc.
USA

Tel: +1 508 880 8990

Bioline GmbH
GERMANY

Tel: +49 (0)337 168 1229

Bioline (Aust) Pty. Ltd
AUSTRALIA

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763

Bioline France
FRANCE

Tel: +33 (0)1 42 56 04 40
Fax: +33 (0)9 70 06 62 10

Meridian Bioscience Asia Pte Ltd
SINGAPORE

Tel: +65 6774 7196
Fax: +65 6774 6441