

PURIFYING TOTAL RNA FROM MICRODISSECTED CRYOSECTIONS, LASER CAPTURED CELLS OR CULTURED CELLS**1 Provide sample**

Transfer sample e.g. microdissected tissue cryosection, pelleted cultured cells (up to 5×10^5), or laser captured cells to a sterile 1.5ml microcentrifuge tube (not supplied).

2 Cell lysis and homogenization

Add 100 μ l Lysis Buffer RLY and 2 μ l TCEP to sample and vortex vigorously (2 x 5s).

3 Add Carrier RNA

Add 5 μ l (20ng) Carrier RNA working solution to lysate.

Mix by vortexing (2 x 5s).

Briefly spin down (1s at 1000 x g).

4 Filter lysate (optional)

Place ISOLATE II Filter (violet) in a 2ml Collection Tube (supplied).

Load lysate and centrifuge 30s at 11,000 x g.

Discard ISOLATE II Filter.

Step 4 may be omitted when processing small amounts of sample, e.g. $<10^6$ cells.

5 Adjust RNA binding conditions

Add 100 μ l ethanol (70%) to homogenized lysate.

Mix by pipetting up and down 5 times.

6 Bind RNA

Place ISOLATE II RNA Micro Column (blue) in a 2ml Collection Tube.

Load lysate onto column and centrifuge 30s at 11,000 x g.

Place column in a new 2ml Collection Tube.

7 Desalt silica membrane

Add 100 μ l Membrane Desalting Buffer (MEM).

Centrifuge 30s at 11,000 x g to dry membrane.

Re-use Collection Tube.

8 Digest DNA

Add 3 μ l reconstituted DNase I to 27 μ l Reaction Buffer for DNase I (RDN).

Mix by gently flicking tube.

Apply 25 μ l DNase I reaction mixture directly onto center of silica membrane.

Incubate at room temperature for 15 min.

9 Wash and dry silica membrane**1st Wash**

- Add 100 μ l Wash Buffer RW1.

Incubate for 2 min at room temperature.

Centrifuge 30s at 11,000 x g.

Place column into a new 2ml Collection Tube.

2nd Wash

- Add 400 μ l Wash Buffer RW2.

Centrifuge 30s at 11,000 x g.

Discard flow-through and place column back into Collection Tube.

3rd Wash

- Add 200 μ l Wash Buffer RW2.

Centrifuge 2 min at 11,000 x g to dry membrane completely.

Place column into a nuclease-free 1.5ml Collection Tube (supplied).

10 Elute RNA

Add 10 μ l RNase-free water (supplied) directly onto center of silica membrane.

Centrifuge at 11,000 x g for 30s.