

## Associated products

| Product                       | Description                                                                                        | Pack Size                         | Cat No.                             |
|-------------------------------|----------------------------------------------------------------------------------------------------|-----------------------------------|-------------------------------------|
| ISOLATE II Genomic DNA Kit    | Rapid isolation of high-quality genomic DNA from many different starting material                  | 10 Preps<br>50 Preps<br>250 Preps | BIO-52065<br>BIO-52066<br>BIO-52067 |
| ISOLATE II Plant DNA Kit      | Rapid isolation of high-quality genomic DNA from a wide variety of plant species                   | 10 Preps<br>50 Preps<br>250 Preps | BIO-52068<br>BIO-52069<br>BIO-52070 |
| ISOLATE II RNA Mini Kit       | Isolation of high-yield and extremely pure total RNA from a variety of samples                     | 10 Preps<br>50 Preps<br>250 Preps | BIO-52071<br>BIO-52072<br>BIO-52073 |
| ISOLATE II RNA Plant Kit      | Isolation of high-yield and extremely pure total RNA from a wide variety of plant species          | 10 Preps<br>50 Preps              | BIO-52076<br>BIO-52077              |
| TRIsure™                      | Quick isolation of high-quality RNA from a variety of sources for subsequent use in cDNA synthesis | 100 mL<br>200 mL                  | BIO-38032<br>BIO-38033              |
| SensiFAST™ cDNA Synthesis Kit | Fully optimized to generate maximum yields of full-length and low abundance cDNA from RNA          | 50 Reactions<br>250 Reactions     | BIO-65053<br>BIO-65054              |
| Agarose                       | Molecular biology grade agarose                                                                    | 100 g<br>500 g                    | BIO-41026<br>BIO-41025              |

## Technical support

If the troubleshooting guide does not solve the difficulty you are experiencing, please contact Technical Support with details of reaction setup, cycling conditions and relevant data.

Email: [mbi.tech@meridianlifescience.com](mailto:mbi.tech@meridianlifescience.com)

## Trademark and licensing information

1) Trademarks: SensiFAST™ (Bioline Reagents Ltd), SYBR® (Molecular Probes), ROX™, StepOne™ (ABI), Mx4000, Mx3000P and Mx3005P (Stratagene), iCycler™, MyiQ5™, Opticon™, Chromo4™, MiniOpticon™, (Bio-Rad), LightCycler®, TaqMan® (Roche), SmartCycler™ (CEPheid), RotorGene™, Scorpion® (Qiagen), RealPlex™ (Eppendorf), Quantica™ (Techne)

Bioline Reagents Ltd  
UNITED KINGDOM  
Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

Bioline USA Inc.  
USA  
Tel: +1 901 382 8716  
Fax: +1 901 382 0027

Bioline GmbH  
GERMANY  
Tel: +49 (0)3371 60222 00  
Fax: +49(0)3371 60222 01

Bioline (Aust) Pty. Ltd  
AUSTRALIA  
Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

## SensiFAST™ Probe No-ROX One-Step Kit

Shipping: On dry/blue ice Catalog numbers

Batch No.: See vial BIO-76001: 100 x 20 µL reactions: 1 x 1 mL

Concentration: see vial BIO-76005: 500 x 20 µL reactions: 5 x 1 mL

**bioline**  
meridian BIOSCIENCE

Store at -20°C

### Storage and stability:

The SensiFAST Probe No-ROX One-Step Kit is shipped on dry/blue ice. All kit components should be stored at -20°C upon receipt. Excessive freeze/thawing is not recommended.

### Expiry:

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

### Quality control:

The SensiFAST Probe No-ROX One-Step Kit and its components are extensively tested for activity, processivity, efficiency, heat activation, sensitivity, absence of nuclease contamination and absence of nucleic acid contamination.

### Safety precautions:

Please refer to the material safety data sheet for further information.

### Notes:

This reagent has been manufactured under 13485 Quality Management System for research and/or further manufacturing use.

## Description

The SensiFAST™ Probe No-ROX One-Step Kit has been formulated for highly reproducible first-strand cDNA synthesis and subsequent real-time PCR in a single tube. The kit is formulated for use with probe-detection technology, including TaqMan®, Scorpions® and molecular beacon probes. A combination of the latest advances in buffer chemistry together with a reverse transcriptase and hot-start DNA polymerase system ensures that SensiFAST Probe No-ROX One-Step Kit produces fast, highly-specific and ultra-sensitive one-step real-time RT-PCR.

The SensiFAST Probe No-ROX One-Step Kit consists of a 2x SensiFAST Probe One-Step mix, separate reverse transcriptase and RiboSafe RNase Inhibitor.

## Kit components

| Reagent                                   | 100 x 20 µL reactions | 500 x 20 µL reactions |
|-------------------------------------------|-----------------------|-----------------------|
| SensiFAST™ Probe No-ROX One-Step mix (2x) | 1 x 1 mL              | 5 x 1 mL              |
| RiboSafe RNase Inhibitor                  | 1 x 40 µL             | 1 x 200 µL            |
| Reverse transcriptase                     | 1 x 20 µL             | 1 x 100 µL            |
| DEPC-H <sub>2</sub> O                     | 1 x 1.8 mL            | 2 x 1.8 mL            |

## Instrument compatibility

The SensiFAST Probe No-ROX One-Step Kit has been optimized for use with all probe chemistries, including TaqMan, FRET, Scorpions and molecular beacon probes.

The SensiFAST Probe No-ROX One-Step Kit can be used on all real-time PCR instruments.

## General considerations

When handling RNA, it is important to use RNase-free plasticware and reagents. We also recommend performing RNA work in an RNase-free area. To help prevent any carry-over DNA contamination, we recommend that separate areas are maintained for reaction set-up, PCR amplification and any post-PCR gel analysis. It is essential that any tubes containing amplified DNA product are not opened in the reaction set-up area.

**Primers and probe:** These guidelines refer to the use of dual-labeled probes. Please refer to the relevant literature when using other probe types. The sequence and concentration of the probe and primers, as well as amplicon length, can be critical for specific amplification, yield and overall efficiency of any real-time RT-PCR.

We strongly recommend taking the following points into consideration when designing and running your real-time RT-PCR:

- use primer-design software, such as Primer3 (<http://frodo.wi.mit.edu/primer3/>) or visual OMP™ (<http://dnasoftware.com/>). Primers should have a melting temperature (T<sub>m</sub>) of approximately 60°C. The T<sub>m</sub> of the probe should be approximately 10°C higher than that of the primers
- optimal amplicon length should be 80-200 bp, and should not exceed 400 bp
- final primer concentration of 400 nM is suitable for most probe reactions. However, to determine the optimal concentration we recommend titrating in the range 0.2-1 µM
- use an equimolar primer concentration
- a final probe concentration of 100 nM is suitable for most applications. We recommend that the final probe concentration is at least 2-fold lower than the primer concentration  
*Note: In multiplex real-time RT-PCR, probe concentrations in excess of 100 nM can result in cross-channel fluorescence*
- where possible, use intron-spanning primers to avoid amplification from genomic DNA

**Template:** It is important that the RNA template is intact and devoid of DNA or contaminating inhibitors of both reverse transcription and PCR. For high purity RNA, we recommend using the Bioline ISOLATE II RNA Mini Kit (BIO-52073). RNA stocks and dilutions should be made in DEPC-treated water to avoid any RNase-mediated degradation.

The recommended amount of template for one-step real-time RT-PCR is dependent upon the type of RNA used:

- **total RNA:** purified total RNA can be used in the range from 1 pg to 1 µg per 20 µL reaction
- **mRNA:** purified mRNA can be used from 0.01 pg per 20 µL reaction

**MgCl<sub>2</sub>:** The MgCl<sub>2</sub> concentration in the 1x reaction mix is 3 mM. In the majority of real-time RT-PCR conditions this is optimal for both the reverse transcriptase and the hot-start DNA polymerase. If necessary, we suggest titrating the MgCl<sub>2</sub> to a maximum of 5 mM.

**RT-PCR controls:** It is important to detect the presence of contaminating DNA that may affect the reliability of the data. Always include a no-RT control reaction, by omitting the reverse transcriptase from the reaction.

## Procedure

**Reaction mix composition:** Prepare a real-time RT-PCR mastermix. The volumes given below are based on a standard 20 µL final reaction mix and can be scaled accordingly.

| Reagent                                | Volume      | Final concentration |
|----------------------------------------|-------------|---------------------|
| 2x SensiFAST Probe No-ROX One-Step Mix | 10 µL       | 1x                  |
| 10 µM Forward Primer                   | 0.8 µL      | 400 nM              |
| 10 µM Reverse Primer                   | 0.8 µL      | 400 nM              |
| 10 µM Probe                            | 0.2 µL      | 100 nM              |
| Reverse transcriptase                  | 0.2 µL      | -                   |
| RiboSafe RNase Inhibitor               | 0.4 µL      | -                   |
| H <sub>2</sub> O                       | up to 16 µL |                     |
| Template                               | 4 µL        |                     |
| <b>20 µL Final volume</b>              |             |                     |

**Suggested RT-qPCR conditions:** The following real-time RT-PCR conditions are suitable for the SensiFAST Probe No-ROX One-Step Kit with the majority of amplicons and real-time PCR instruments. However, the cycling conditions can be varied to suit different probe-based reactions or machine-specific protocols. The detection channel on the real-time instrument should be set to acquire at the appropriate wavelength(s). We recommend using the following cycling conditions for optimal results:

**Sensitivity testing and C<sub>t</sub> values:** When comparing SensiFAST with a mix from another supplier we strongly recommend amplifying from a 10-fold template dilution series. Loss of detection at low template concentration is the only direct measurement of sensitivity. An early C<sub>t</sub> value is not an indication of good sensitivity, but rather an indication of speed.

### • Cycling for dual-labeled probes

| Cycles | Temp.        | Time        | Notes                                                        |
|--------|--------------|-------------|--------------------------------------------------------------|
| 1      | 45°C         | 10 min      | Reverse transcription                                        |
| 1      | 95°C         | 2 min       | Polymerase activation                                        |
| 40     | 95°C<br>60°C | 5 s<br>20 s | Denaturation<br>Annealing/extension (acquire at end of step) |

**Real-Time RT-PCR optimization:** The following optimization may be necessary to improve the efficiency of some reactions, such as multiplexing with more than two probes, or if the target amplicon is longer than 200 bp.

- The reverse transcription reaction time can be extended up to 20 minutes and/or the temperature can be increased up to 48°C
- The annealing/extension time can be extended up to 60 seconds and/or the temperature can be increased up to 65°C

## Troubleshooting guide (Continued)

| Problem                                                    | Possible Cause                                              | Recommendation                                                                                                                                         |
|------------------------------------------------------------|-------------------------------------------------------------|--------------------------------------------------------------------------------------------------------------------------------------------------------|
| Non-specific amplification product<br>AND<br>Primer-dimers | Inefficient reverse transcription                           | Extend reverse transcription time up to 20 min and/or increase the temperature up to 48°C                                                              |
|                                                            | Suboptimal primer/probe design                              | Redesign primers and/or probe using appropriate software, or use validated primers/probes                                                              |
|                                                            | Primer/probe concentration too high                         | Test dilution series of primer/probe concentrations until primer-dimer/non-specific amplification products disappear                                   |
|                                                            | Primer/probe concentration too low                          | Use primer concentration between 200 nM and 1 µM and probe concentration at least 2 fold lower                                                         |
|                                                            | Annealing/extension temperature too low                     | Increase annealing/extension temperature up to 65°C or until primer-dimer/non-specific amplification products disappear                                |
|                                                            | Template concentration too low                              | Increase template concentration                                                                                                                        |
|                                                            | Template concentration too high                             | Reduce template concentration until non-specific products disappear                                                                                    |
|                                                            | Extension time too long                                     | Reduce extension time to determine whether non-specific products are reduced                                                                           |
| Variability between replicates                             | Error in reaction set-up                                    | Prepare large volume mastermix                                                                                                                         |
|                                                            | Air bubbles in reaction mix                                 | Centrifuge reaction samples/plate prior to running on a real-time instrument                                                                           |
| Late amplification trace                                   | Inefficient reverse transcription                           | Extend reverse transcription time up to 20 min and/or increase the temperature up to 48°C                                                              |
|                                                            | Activation time too short                                   | Ensure SensiFAST Probe No-ROX One-Step mix is activated for a minimum of 1 min at 95°C before cycling                                                  |
|                                                            | Annealing temperature too high                              | Decrease annealing temperature in steps of 2°C                                                                                                         |
|                                                            | Extension time too short                                    | Double extension time to determine whether the cycle threshold (C <sub>T</sub> ) is affected                                                           |
|                                                            | Template concentration too low                              | Increase concentration if possible                                                                                                                     |
|                                                            | Template is degraded                                        | Re-isolate template from sample material or use freshly prepared template dilution                                                                     |
|                                                            | Suboptimal primer/probe design                              | Redesign primers/probe using appropriate software, or use validated primers                                                                            |
|                                                            | Primer/probe concentration too low                          | Increase concentration of primers in 100 nM increments and probe concentration in increments at least 2 fold lower than that of the primer             |
| RNase contamination                                        | Ensure RNase inhibitor is added before addition of template |                                                                                                                                                        |
| PCR efficiency below 90%                                   | Extension time too short                                    | Increase extension time                                                                                                                                |
|                                                            | Primer concentration too low                                | Increase concentration of primers in 100nM increments                                                                                                  |
|                                                            | Suboptimal primer/probe design                              | Redesign primer/probe using appropriate software or use validated primer/probe                                                                         |
| PCR efficiency above 110%                                  | Template is degraded or contains PCR inhibitors             | Re-isolate template from sample material, or use freshly prepared template dilution, or purify template and resuspend it in water                      |
|                                                            | Non specific amplification and/or primer-dimers             | Use 4% agarose gel electrophoresis to confirm presence of non-specific amplification products. See above for preventing/removing non-specific products |

## Troubleshooting guide

| Problem                                                             | Possible Cause                                                                         | Recommendation                                                                                                                                                                                                                                                                                     |
|---------------------------------------------------------------------|----------------------------------------------------------------------------------------|----------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| No amplification trace<br>AND<br>No product on agarose gel          | Activation time too short                                                              | Ensure SensiFAST Probe No-ROX One-Step mix is activated for a minimum of 2 min at 95°C before cycling                                                                                                                                                                                              |
|                                                                     | Error in protocol setup                                                                | Verify that correct reagent concentrations, volumes, dilutions and storage conditions have been used                                                                                                                                                                                               |
|                                                                     | Suboptimal primer/probe design                                                         | Use primer design software or validated primers/probes. Test assay on a control template                                                                                                                                                                                                           |
|                                                                     | Incorrect concentration of primers/probe                                               | Use primer concentrations between 200 nM and 1 µM. Probe concentration should be at least 2-fold lower than the primer concentration                                                                                                                                                               |
|                                                                     | Template degraded                                                                      | Re-isolate your template from the sample material or use freshly prepared template dilution. We recommend using the ISOLATE RNA kits for template preparation and DEPC-treated water for resuspension or dilution of the template<br>Verify the integrity of RNA using agarose gel electrophoresis |
|                                                                     | Primers/probe degraded                                                                 | Use newly synthesized primers and/or probe                                                                                                                                                                                                                                                         |
|                                                                     | Template contaminated with real-time RT-PCR inhibitors                                 | Further dilute template before real-time RT-PCR or purify template and resuspend it in DEPC-treated water                                                                                                                                                                                          |
|                                                                     | Template concentration too low                                                         | Increase concentration used                                                                                                                                                                                                                                                                        |
| Cycling conditions not optimal                                      | Increase extension/annealing time, increase cycle number, reduce annealing temperature |                                                                                                                                                                                                                                                                                                    |
| No amplification trace<br>AND<br>PCR product present on agarose gel | Error in instrument setup                                                              | Check that the acquisition settings are correct during cycling                                                                                                                                                                                                                                     |