

## **PURIFYING GENOMIC DNA FROM CULTURED CELLS AND HUMAN OR ANIMAL TISSUE**

---

### **1 SAMPLE PREPARATION**

#### **1.1 HUMAN OR ANIMAL TISSUE**

Cut up 25 mg tissue and transfer to 1.5 mL microcentrifuge tube (proceed to step 2).

#### **1.2 CULTURED CELLS**

Resuspend up to  $10^7$  cells in 200  $\mu$ L Lysis Buffer GL.

Add 25  $\mu$ L Proteinase K solution and 200  $\mu$ L Lysis Buffer G3.

Incubate at 70°C for 10-15 min (proceed to step 4).

### **2 PRE-LYSIS**

Add 180  $\mu$ L Lysis Buffer GL and 25  $\mu$ L Proteinase K solution.

Completely cover sample with solution and vortex.

Incubate at 56°C for 1–3 hours (until completely lysed), shake or vortex occasionally.

### **3 LYSE SAMPLE**

Vortex sample briefly and add 200  $\mu$ L Lysis Buffer G3.

Vortex vigorously and incubate at 70°C for 10 min.

### **4 ADJUST DNA BINDING CONDITIONS**

Vortex briefly and add 210  $\mu$ L ethanol (96-100%) to sample.

Vortex vigorously.

### **5 BIND DNA**

Place ISOLATE II Genomic DNA Spin Column (green) in a 2 mL Collection Tube.

Load sample to column and centrifuge 1 min at 11,000 x g.

Discard flow-through and reuse Collection Tube.

### **6 WASH SILICA MEMBRANE**

- Add 500  $\mu$ L Wash Buffer GW1.

Centrifuge 1 min at 11,000 x g.

Discard flow-through and reuse Collection Tube.

- Add 600  $\mu$ L Wash Buffer GW2.

Centrifuge 1 min at 11,000 x g.

Discard flow-through and reuse Collection Tube.

### **7 DRY SILICA MEMBRANE**

Centrifuge 1 min at 11,000 x g, to remove residual ethanol.

Place ISOLATE II Genomic DNA Spin Column in a 1.5 mL microcentrifuge tube (not supplied).

### **8 ELUTE DNA**

Add 100  $\mu$ L preheated Elution Buffer G (70°C) onto center of silica membrane.

Incubate at room temperature for 1 min.

Centrifuge 1 min at 11,000 x g.