

PURIFYING TOTAL RNA FROM CULTURED CELLS AND TISSUE

1 SAMPLE HOMOGENIZATION

See manual for recommended sample amounts and homogenization methods.

2 CELL LYSIS

Add 350 μ L Lysis Buffer RLY and 3.5 μ L β -ME to cell pellet (up to 5×10^6 eukaryotic cultured cells) or to ground tissue (up to 30 mg tissue) and vortex vigorously.

3 FILTER LYSATE

Place ISOLATE II Filter (violet) in a 2 mL Collection Tube (supplied).

Load lysate and centrifuge 1 min at 11,000 x g.

If a visible pellet forms, transfer supernatant avoiding any pellet to a new 1.5 mL microcentrifuge tube (not supplied).

Alternatively, pass lysate 5-10 times through a nuclease-free 20 gauge (0.9 mm) needle and syringe.

4 ADJUST RNA BINDING CONDITIONS

Discard ISOLATE II Filter and add 350 μ L ethanol (70%) to homogenized lysate.

Mix by pipetting up and down 5 times.

Alternatively, transfer flow-through into a new 1.5 mL microcentrifuge tube (not supplied), add 350 μ L ethanol (70%) and mix by vortexing (2 x 5s).

5 BIND RNA

Place ISOLATE II RNA Mini Column (blue) in a 2 mL Collection Tube.

Pipette lysate up and down 2–3 times and load lysate onto column.

Centrifuge 30s at 11,000 x g.

Place column in a new 2 mL Collection Tube.

6 DESALT SILICA MEMBRANE

Add 350 μ L Membrane Desalting Buffer (MEM).

Centrifuge at 11,000 x g for 1 min to dry membrane.

7 DIGEST DNA

Add 10 μ L reconstituted DNase I to 90 μ L Reaction Buffer for DNase I (RDN).

Mix by gently flicking tube.

Apply 95 μ L DNase I reaction mixture directly onto center of silica membrane.

Incubate at room temperature for 15 min.

8 WASH AND DRY SILICA MEMBRANE

1st Wash

- Add 200 μ L Wash Buffer RW1.

Centrifuge 30s at 11,000 x g.

Place column into a new 2 mL Collection Tube.

2nd Wash

- Add 600 μ L Wash Buffer RW2.

Centrifuge 30s at 11,000 x g.

Discard flow-through and place column back into Collection Tube.

3rd Wash

- Add 250 μ L Wash Buffer RW2.

Centrifuge 2 min at 11,000 x g to dry membrane completely.

Place column into a nuclease-free 1.5 mL Collection Tube (supplied).

9 ELUTE RNA

Add 60 μ L RNase-free water (supplied) directly onto center of silica membrane.

Centrifuge at 11,000 x g for 1 min.