

PURIFYING TOTAL RNA FROM CULTURED CELLS AND TISSUE**1 Sample homogenization**

See manual for recommended sample amounts and homogenization methods.

2 Cell lysis

Add 350µl Lysis Buffer RLY and 3.5µl β-ME to cell pellet (up to 5x10⁶ eukaryotic cultured cells) or to ground tissue (up to 30mg tissue) and vortex vigorously.

3 Filter lysate

Place ISOLATE II Filter (violet) in a 2ml Collection Tube (supplied).

Load lysate and centrifuge 1 min at 11,000 x g.

If a visible pellet forms, transfer supernatant avoiding any pellet to a new 1.5ml microcentrifuge tube (not supplied).

Alternatively, pass lysate 5-10 times through a nuclease-free 20 gauge (0.9mm) needle and syringe.

4 Adjust RNA binding conditions

Discard ISOLATE II Filter and add 350µl ethanol (70%) to homogenized lysate. Mix by pipetting up and down 5 times.

Alternatively, transfer flow-through into a new 1.5ml microcentrifuge tube (not supplied), add 350µl ethanol (70%) and mix by vortexing (2 x 5s).

5 Bind RNA

Place ISOLATE II RNA Mini Column (blue) in a 2ml Collection Tube.

Pipette lysate up and down 2–3 times and load lysate onto column.

Centrifuge 30s at 11,000 x g.

Place column in a new 2ml Collection Tube.

6 Desalt silica membrane

Add 350µl Membrane Desalting Buffer (MEM).

Centrifuge at 11,000 x g for 1 min to dry membrane.

7 Digest DNA

Add 10µl reconstituted DNase I to 90µl Reaction Buffer for DNase I (RDN).

Mix by gently flicking tube.

Apply 95µl DNase I reaction mixture directly onto center of silica membrane.

Incubate at room temperature for 15 min.

8 Wash and dry silica membrane

1st Wash

- Add 200µl Wash Buffer RW1.
Centrifuge 30s at 11,000 x g.
Place column into a new 2ml Collection Tube.

2nd Wash

- Add 600µl Wash Buffer RW2.
Centrifuge 30s at 11,000 x g.
Discard flow-through and place column back into Collection Tube.

3rd Wash

- Add 250µl Wash Buffer RW2.
Centrifuge 2 min at 11,000 x g to dry membrane completely.
Place column into a nuclease-free 1.5ml Collection Tube (supplied).

9 Elute RNA

Add 60µl RNase-free water (supplied) directly onto center of silica membrane.

Centrifuge at 11,000 x g for 1 min.