

PURIFYING GENOMIC DNA FROM WHOLE BLOOD**1 Lyse blood**

Add 25µl Proteinase K solution and 200µl sample into a 1.5ml microcentrifuge tube (not supplied).

Add 200µl Lysis Buffer G3 and vortex vigorously for 10–20s.

Incubate samples at 70°C for 10–15 min.

2 Adjust DNA binding conditions

Add 210µl ethanol (96–100%) and vortex.

3 Bind DNA

Place ISOLATE II Blood DNA Spin Column in a 2ml Collection Tube and load sample onto column.

Centrifuge 1 min at 11,000 x g.

Repeat at a higher g force if samples are not completely filtered through matrix.

Place column in a new 2ml Collection Tube.

4 Wash silica membrane

- Add 500µl Wash Buffer GW1.

Centrifuge 1 min at 11,000 x g.

Place column in a new 2ml Collection Tube.

- Add 600µl Wash Buffer GW2.

Centrifuge 1 min at 11,000 x g.

Discard flow-through and reuse Collection Tube.

5 Dry silica membrane

Centrifuge 1 min at 11,000 x g to remove residual ethanol.

Place ISOLATE II Blood DNA Spin Column in a 1.5ml microcentrifuge tube (not supplied).

6 Elute DNA

Add 100µl preheated Elution Buffer G (70°C) directly onto center of silica membrane.

Incubate at room temperature for 1 min.

Centrifuge 1 min at 11,000 x g.