5x MyTaq™ Reaction Buffer Colorless

Shipping: On Dry/Blue Ice

Catalog numbers

Batch No.: See vial BIO-37111: 4 x 1 mL

Store at -20 °C



Storage and stability: 5x MyTaq Reaction Buffer Colorless is shipped on Dry/Blue Ice and stored at $-20\,^{\circ}\text{C}$.

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

Safety precautions:

This product is for R&D use only, not for human use, or any other use. Irritating to eyes, respiratory system and skin. May cause sensitization by inhalation and skin contact. Please refer to the material safety data sheet for information regarding hazards and safe handling practice.

Quality control specifications:

5x MyTaq Reaction Buffer Colorless is extensively tested for the absence of nucleic acid contamination prior to release.

Notes:

For research or further manufacturing use only.

Trademarks:

MyTaq is a trademark of Bioline Reagents Ltd.

Description

5x MyTaq™ Reaction Buffer Colorless is an advanced formulation buffer that saves time and delivers superior results, containing dNTPs, MgCl₂ and enhancers at optimal concentrations which eliminates the need for optimization.

Applications

- For use in reactions containing MyTaq DNA Polymerase
- Use at 1 x concentration in reaction mix

Components

Product Name	Volume
5x MyTaq Reaction Buffer Colorless	4 x 1 mL

5x MyTag Reaction Buffer Colorless: The buffer comprises 5 mM dNTPs, 15 mM MgCl₂, stabilizers and enhancers. The concentration of each component has been extensively optimized, reducing the need for further optimization. Additional PCR enhancers such as dextran sulphate or betaine etc. are not recommended.

Protocol

Add 10 µL to a 50 µL PCR reaction

Associated products

Product Name	Pack size	Cat. No.
MyTaq™ DNA polymerase	500 Units	BIO-21105
MyTaq™ HS DNA polymerase	250 Units	BIO-21111

Bioline Ltd UNITED KINGDOM

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822 Meridian Life Science Inc.

Tel: +1 901.382.8716

Bioline GmbH

Tel: +49 (0)3371 60222 03 Fax: +49 (0)3371 60222 01 Bioline (Aust) Pty. Ltd

Tel: +61 (0)2 9209 4180 Fax: +61 (0)2 9209 4763