

Certificate of Analysis

COA No: CA_XBE-0001

Version: 08

BIOTAQ™ DNA Polymerase

For Research and Further Manufacturing use only

| Catalog No: | BIO-21060 |
|---------------------|---------------|
| Lot No: | PL313-B124090 |
| Storage Conditions: | -20°C |
| Component Lot No: | BT-224101B |
| Expiry date: | February 2026 |

Quality Control Parameters

| Analysis | Specification | Result |
|--|--|--------|
| Functional | A 3Kb fragment is amplified with a dilution series of human genomic DNA and a dilution series of enzyme, using standard conditions and 30 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |
| DNA contamination | Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample. | Passed |
| DNase contamination Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase. | | Passed |

QA / QC Representative:

J. Rahnenführer

Date: 7th February 2024

United Kingdom

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822

USA Tel: +1 901.382.8716 <u>Germany</u>

Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01



Certificate of Analysis

COA No: CA_XBB-0014

Version: 09

MgCl₂ Solution, 50mM

For research or further manufacturing use only

| Catalog No: | BIO-21060 |
|---------------------|---------------|
| Lot No: | PL313-B124090 |
| Storage Conditions: | -20°C |
| Component Lot No: | MG-2031.018 |
| Expiry date: | February 2026 |

Quality Control Parameters

| Analysis | Specification | Result |
|--|---|--------|
| Functional | Fragments of sizes 800bp and 3000bp are amplified with a dilution series of BIOTAQ™ DNA Polymerase, using standard conditions and 30 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |
| DNA contamination | Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample. | Passed |
| DNase contamination Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10 ⁻³ U DNase. | | Passed |

QA / QC Representative:

J. Rahnenführer

Date: 7th February 2024

United Kingdom Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822

<u>USA</u> Tel: +1 901.382.8716 Fax: +1 901.382.0027

Germany

Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01



Certificate of Analysis

COA No: CA XBB-0001

Version: 08

NH₄ Buffer, 10x

For research or further manufacturing use only

| Catalog No: | BIO-21060 | |
|---------------------|---------------|--|
| Lot No: | PL313-B124090 | |
| Storage Conditions: | -20°C | |
| Component Lot No: | NH-324101B | |
| Expiry date: | February 2026 | |

Quality Control Parameters

| Analysis | Specification | Result |
|---|---|--------|
| Functional | Fragments of sizes 800bp and 3000bp are amplified with a dilution series of BIOTAQ TM DNA Polymerase, using standard conditions and 30 cycles. Fragment of size 3000bp is amplified with a dilution series of Lambda DNA template, using standard conditions and 30 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained). | Passed |
| DNA contamination | Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with a reference sample. | Passed |
| Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5 x 10-3 U DNase. | | Passed |

QA / QC Representative:

J. Rahnenführer

Date: 7th February 2024

United Kingdom

Tel: +44 (0)20 8830 5300 Fax: +44 (0)20 8452 2822

<u>USA</u> Tel: +1 901.382.8716 Germany

Tel: +49 (0)3371 60222 00 Fax: +49 (0)3371 60222 01