

	Certificate of Analysis	COA No: CA_RES-0077
		Version: 04

JetSeq Clean For Research Use Only	Catalog No:	BIO-68030
	Lot No:	NG874-B080330
	Shipping / Storage Conditions:	at +4°C
	Component Lot No:	NG874-B080330
	Expiry date:	March 2022

Quality Control Parameters

Analysis	Specification	Result
Small fragment removal test	JetSeq Clean beads are able to remove 50 bp fragments from a population of DNA fragments of various sizes, in standard clean-up conditions.	PASS
Double-sided size selection	JetSeq Clean beads are able to select a specific size range from a mixture of DNA sizes, in standard double-sided size selection conditions.	PASS

Authorised by Christopher Weatherall



United Kingdom
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

info@meridianlifescience.com
Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

info.de@bioline.com
Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

France

info.fr@bioline.com
Tel: +33 (0)1 42 56 04 40
Fax: +33 (0)9 70 06 62 10

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763