

## Components

Tetro cDNA Synthesis Kit	30 Reactions	100 Reactions
Cat No.	BIO-65042	BIO-65043
5x RT Buffer	120µl	400µl
(200u/µl) Reverse Transcriptase	30µl	100µl
(10u/µl) RiboSafe RNase Inhibitor	30µl	100µl
dNTP Mix 10mM Total	30µl	100µl
Oligo (dT) <sub>18</sub> Primer Mix	30µl	100µl
Random Hexamer Primer Mix	30µl	100µl
DEPC-treated Water	1.8ml	1.8ml

[www.bioline.com](http://www.bioline.com)

