

5x MangoTaq™ Colored Reaction Buffer

Shipping: On Dry/Blue Ice Catalog numbers
Exp. Date: See vial BIO-37102 : 10ml
Batch No.: See vial
Concentration: 5x



A Meridian Life Science® Company

Store at -20°C

Storage and stability:

5x MangoTaq Colored Buffer is shipped on Dry/Blue Ice and can be stored for up to 12 months at -20°C.

Safety precautions:

Harmful if swallowed. Irritating to eyes, respiratory system and skin. Please refer to the material safety data sheet for further information.

Notes:

MangoTaq, MangoMix and HyperLadder are Trademarks of Bioline Limited
This product insert is a declaration of analysis at the time of manufacture.
Research Use Only.

Features

- Optimized reaction buffer for PCR
- Two inert dyes for easy orientation

Applications

- To monitor migration rate during agarose electrophoresis
- For high-throughput application in combination with MangoTaq

Description

The 5x MangoTaq™ Colored Reaction Buffer contains two inert dyes (red and orange) that separate during gel electrophoresis and provide quick reference points for monitoring the mobility of the DNA samples in the gel. The red dye migrates at 300bp on a 2% agarose gel. The orange dye migrates at 30bp on a 2% agarose gel and <10bp on a 1% gel.

The colored reaction buffer can be loaded directly onto an agarose gel for analysis, without the need for separate gel-loading buffer. The presence of the dyes has no effect on routine enzymatic manipulations, although rare exceptions may exist. For high-throughput applications, the colored reaction buffer and MangoTaq DNA Polymerase make an ideal choice, since this combination enables the user to load the PCR products directly on a gel to facilitate visual recognition

Components

Product Name	Volume
5x MangoTaq™ Colored Reaction Buffer	10 ml

Associated products

Product Name	Pack size	Cat. No.
MangoTaq™	1000 Units	BIO-21083
MangoTaq™	2000 Units	BIO-21082
MangoMix™	250 Reactions	BIO-25033
MangoMix™	1000 Reactions	BIO-25034
HyperLadder™ I	200 Lanes	BIO-33025
Agarose	100g	BIO-41026
SureClean Plus	1 x 5ml	BIO-37047

Product Citations:

1. Gillespie, K., Valovin, S., Saunby, J., Hunter, K., Savage, D., Middleton, D., Todd, J., Bingley, P., Gale, E.A.M. *Tissue Antigens*, **56(6)**, 530 (2000)
2. Caldwell, G., Jones, C., Gensberg, K., Jan, S., Hardy, R., Byrd, P., Chughtai, S., Wallis, Y., Matthews, G. and Morton, D.G. *Cancer Research*, **64**, 883-888 (2004)
3. Al-Sharif, F., Makki, R., Ollier, W., Hajeer, A. *Journal Immunogenetics*, **49(7-8)**, 740-741 (1999)

Bioline Ltd
UNITED KINGDOM

Tel: +44(0)20 8830 5300
Fax: +44 (0)20 8452 2822

Bioline USA Inc.
USA

Tel: +1 508 880 8990
Fax: +1 508 880 8993

Bioline GmbH
GERMANY

Tel: +49(0)33 7168 1229
Fax: +49 (0)337168 1244

Bioline (Aust) Pty. Ltd
AUSTRALIA

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763