



Certificate of Analysis

COA No: CA_RES-0019

Version: 01

ISOLATE II Plasmid Mini Kit

For Research Use Only

Storage Conditions: Room Temperature

Lot number: IS494:314I

Expiry date: N/A

Quality Control Parameters

| Component | Specification | Result | Batch |
|-----------------------------------|---|--------|----------|
| ISOLATE II Plant DNA Spin Columns | Nucleic acid recovery >90% | Passed | SC-314I |
| Collection tubes | n/a | n/a | CT-314I |
| Resuspension Buffer P1 | pH = 7.9-8.1, density | Passed | P1-314I |
| Lysis Buffer P2 | density | Passed | P2-314I |
| Neutralization Buffer P3 | pH = 4.1-4.4, density | Passed | P3-314I |
| Wash Buffer PAW1 | pH = 7.2-7.6, density | Passed | PW1-314I |
| Wash Buffer PAW2 | pH = 7.5-7.9, density, sterility | Passed | PW2-314I |
| Elution Buffer P | pH = 8.3-8.7, density, sterility | Passed | EBP-314I |
| RNase A | No DNase activity detectable in connection with our kit | Passed | RN-314I |
| | | | |
| | | | |

Authorised by Jade James

Europe
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Europe
Germany

info.de@bioline.com
Tel: +49 (0)3371 681 229
Fax: +49 (0)3371 681 244

America

info.us@bioline.com
Tel: +1 508 880 8990
Fax: +1 508 880 8993

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763