



Product Insert
KCl Reaction Buffer

Product:
10x KCl Reaction Buffer

Catalogue Number:
BIO-37066 3 x 1.2ml

Application

- For use in reactions containing DNA Polymerase

Recommendations for Use:
Use at 1x concentration in reaction mix. This will give a $MgCl_2$ concentration of 1.5mM.

Product Specifications

Storage Conditions:
10x KCl Reaction Buffer can be stored for 12 months at -20°C. Avoid multiple freeze/thaw cycles.

Shipping Conditions:
On Dry Ice or Blue Ice

Composition (10x):
500mM KCl, 100mM Tris-HCl (pH 8.8 at 25°C), 15mM $MgCl_2$ and stabilizer.

Associated Products:

Product Name	Pack Size	Cat No
ACCUZYME DNA Polymerase	500 Units	BIO-21052
AccuSure DNA Polymerase	500 Units	BIO-21069
IMMOLASE DNA Polymerase	250 Units	BIO-21046
Diamond DNA Polymerase	250 Units	BIO-21058

Product Citations:

- Santos, F.R., Pandey, A.A., Tyler-Smith, C., Pena, S.D.J., Schanfield, M., Leonard, W.R., Osipoval, L.; Crawford, M.H.; Mitchell, R.J. *Am. j. hum. genet.*, **64(2)**, 619-628, 1999
- Zito, I., Thiselton, D.L., Gorin, M.B., Stout, J.T., Plant, C., Bird, A.C., Bhattacharya, S.S., Hardcastle, A.J. *Journal Human Genetics*, **105(1-2)**, 57-62, 1999
- Bensasson, D., Zhang D-X. and Hewitt, G.M. *Molecular Biology and Evolution* **17**, 406-415, 2000.

Notes

- This product insert is a declaration of analysis at the time of manufacture.
- Research Use Only.

UK
Bioline
16 The Edge Business Centre
Humber Road
London, NW2 6EW
U.K

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Germany
Bioline GmbH
Im Biotechnologiepark TGZ 2
D-14943 Luckenwalde
Germany

Tel: +49 (0)33 7168 1229
Fax: +49 (0)33 7168 1244

email: info@bioline.com
website: www.bioline.com

USA
Bioline USA Inc.
305 Constitution Dr.
Taunton
MA 02780
USA

Toll Free: 888 257 5155
Tel: 508 880 8990
Fax: 508 880 8993

Australia
Bioline (Aust) Pty Ltd
PO Box 122
Alexandria NSW 1435
Australia

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763