

# IPTG (Dioxane-Free)

Shipping: On Dry/Blue Ice      Catalog numbers

Batch No.: See vial      BIO-37036 : 5 g



Store at -20°C

## Storage and stability:

The IPTG is shipped on dry/blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided.

## Expiry

When stored under the recommended conditions and handled correctly, full activity of the kit is retained until the expiry date on the outer box label.

**Safety precautions:** Harmful if swallowed. Irritating to eyes and skin. Please refer to the material safety data sheet for further information and safety handling practice.



Xi: IRRITANT

## Quality control specifications:

IPTG is extensively tested for activity and absence of contamination prior to release.

## Notes:

Research use only.

## Description

Isopropyl- $\beta$ -D-thiogalactopyranoside (IPTG) is a chemical analogue of galactose, which cannot be hydrolysed by the enzyme  $\beta$ -Galactosidase. Hence, it induces the *E. coli lac* operon activity by binding and inhibiting the lac repressor without being degraded. Genes controlled by the *lac* or *tac* promoter/operator sequences are expressed to high levels in the presence of IPTG.

It is recommended for blue/white screening by *lacZ*  $\alpha$ -complementation with appropriate vectors and host strains, and for protein expression with BL21 (DE3) cells and other lac promoter-derived expression systems.

## Features

- Induces *E. coli lac* operon activity
- >99.6% purity by HPLC
- Available as powder and stabilized stock solution

## Applications

- Blue/white color screening
- Induction of *lac* operon for protein expression
- Genes controlled by the *lac* or *tac* promoter/operator sequences are expressed to high levels in the presence of IPTG

## Characteristics

Empirical formula	C <sub>9</sub> H <sub>18</sub> O <sub>5</sub> S
Molecular weight	238.3
Form	1M Solution in water
Melting point	112-113 °C
Purity	≥99% by HPLC

## Associated Products

Product Name	Pack Size	Catalog No.
X-GAL Powder	1 g	BIO-37035
SureClean Plus	5 mL	BIO-37047

## Recommendations for use

- Dissolve to desired concentration in Water, 18.2M $\Omega$  (DNase/RNase-free).
- Filter sterilize through a 0.22 micron disposable filter.
- Store solution at -20 °C.
- For blue/white screening, 0.1mM final IPTG concentration in LB (Luria Broth) media is recommended

## Product Citations:

1. Ali, F.R., *et al. Meth. Mol. Biol.* **628**, 195-214 (2010).
2. Prabhakar, V., *et al. FEBS Lett.* **583(6)**, 983-991 (2009).
3. Chan, C-H., *et al. Conserv. Gen.* **9(4)**, 1067-1070 (2008).
4. Hamblin, K., *et al. Mol. Microbiol.* **68(6)**, 1395-1405 (2008).
5. Maruta, F., *et al. J. Drug Targeting* **15(4)**, 311-319 (2007).
6. Ross, P.J., *et al. Infect. Immun.* **72(3)**, 1568-1579 (2004).
7. Meng, G. & Fütterer, K. *Nat. Struct. Biol.* **10**, 935-941 (2003).

Bioline Reagents Ltd  
UNITED KINGDOM  
Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

Bioline USA Inc.  
USA  
Tel: +1 800 327 6299  
Fax: +1 901 382 0027

Bioline GmbH  
GERMANY  
Tel: +49 (0)337 168 1229  
Fax: +49 (0)337 168 1244

Bioline (Aust) Pty. Ltd  
AUSTRALIA  
Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763

Bioline France  
FRANCE  
Tel: +33 (0)1 42 56 04 40  
Fax: +33 (0)9 70 06 62 10