



Certificate of Analysis

COA No: CA_BS-0019

Version: 01

Human Genomic DNA

For Research Use Only

Storage
Conditions:

-20°C

Lot number:

HG:313108

Expiry date:

August 2015

Quality Control Parameters

Analysis	Specification	Result
Functionality	23Kb sized fragment is amplified with a dilution series of human genomic DNA, using standard conditions and 35 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
Concentration	Concentration measured spectrophotometrically = 200ng/μl A260/A280 = >1.8<2.0	Passed

Authorised by Jade James

Europe
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Europe
Germany

info.de@bioline.com
Tel: +49 (0)3371 681 229
Fax: +49 (0)3371 681 244

America

info.us@bioline.com
Tel: +1 508 880 8990
Fax: +1 508 880 8993

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763