



## Certificate of Analysis

COA No: CA BS-0004

Version: 01

### Proteinase K Solution

For Research Use Only

Storage Conditions: -20°C

Lot number: PKS:313106

Expiry date: June 2015

### Quality Control Parameters

Analysis	Specification	Result
Activity	Kinetic colorimetric assay using a substrate to establish kinetics when compared to reference standard.	Passed

Authorised by Jade James

**Europe**  
Headquarters UK

[info.uk@bioline.com](mailto:info.uk@bioline.com)  
Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**Europe**  
Germany

[info.de@bioline.com](mailto:info.de@bioline.com)  
Tel: +49 (0)3371 681 229  
Fax: +49 (0)3371 681 244

**America**

[info.us@bioline.com](mailto:info.us@bioline.com)  
Tel: +1 508 880 8990  
Fax: +1 508 880 8993

**Australia**

[info.aust@bioline.com](mailto:info.aust@bioline.com)  
Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763