



Certificate of Analysis

COA No: CA BS-0004

Version: 01

Proteinase K Solution

For Research Use Only

Storage Conditions: -20°C

Lot number: PKS-313110

Expiry date: October 2015

Quality Control Parameters

| Analysis | Specification | Result |
|----------|---|--------|
| Activity | Kinetic colorimetric assay using a substrate to establish kinetics when compared to reference standard. | Passed |

Authorised by Jade James

Europe
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Europe
Germany

info.de@bioline.com
Tel: +49 (0)3371 681 229
Fax: +49 (0)3371 681 244

America

info.us@bioline.com
Tel: +1 508 880 8990
Fax: +1 508 880 8993

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763