



## Certificate of Analysis

COA No: CA\_BS-0019

Version: 01

### Human Genomic DNA

For Research Use Only

Storage Conditions:	-20°C
Lot number:	HG:314103
Expiry date:	March 2016

### Quality Control Parameters

Analysis	Specification	Result
Functionality	23Kb sized fragment is amplified with a dilution series of human genomic DNA, using standard conditions and 35 cycles. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
Concentration	Concentration measured spectrophotometrically = 200ng/μl A260/A280 = >1.8<2.0	Passed

Authorised by Jade James

**Europe**  
Headquarters UK

[info.uk@bioline.com](mailto:info.uk@bioline.com)  
Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**Europe**  
Germany

[info.de@bioline.com](mailto:info.de@bioline.com)  
Tel: +49 (0)3371 681 229  
Fax: +49 (0)3371 681 244

**America**

[info.us@bioline.com](mailto:info.us@bioline.com)  
Tel: +1 508 880 8990  
Fax: +1 508 880 8993

**Australia**

[info.aust@bioline.com](mailto:info.aust@bioline.com)  
Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763