



Certificate of Analysis

COA No: CA_RES-0034

Version: 01

ISOLATE II RNA Mini Kit

For Research Use Only

Storage Conditions:	Room Temperature
Lot number:	IS510-BO28220
Expiry date:	12 months from date of purchase

Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II Filters	n/a	n/a	FIL 115H
ISOLATE II RNA Mini Columns	Nucleic acid recovery >90%	Passed	SC 115H
Collection tubes (2ml)	n/a	n/a	CT2 115H
Collection tubes (1.5ml)	n/a	n/a	CT1.5 115H
Lysis Buffer RLY	pH = 6.7-7.3, density	Passed	RLY 115H
Wash Buffer RW1	pH = 6.9-7.5, density	Passed	RW1 115H
Wash Buffer RW2	pH = 7.5-7.9, density, sterility	Passed	RW2 115H
Membrane Desalting Buffer MEM	pH = 7.8-7.1, density	Passed	MEM 115H
Reaction Buffer for DNase I RDN	pH = 6.9-7.1, density	Passed	RDN 115H
DNase, RNase-free	No RNase activity detectable in connection with our kit	Passed	DN1 115H
RNase-free Water	Density, no RNase activity, sterility	Passed	H2O 115H

Authorised by Jade James

Europe
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Europe
Germany

info.de@bioline.com
Tel: +49 (0)3371 681 229
Fax: +49 (0)3371 681 244

America

info.us@bioline.com
Tel: +1 508 880 8990
Fax: +1 508 880 8993

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763