



Certificate of Analysis

COA No: CA RES-0004

Version: 02

Proteinase K

For Research Use Only

Storage Conditions: -20°C

Lot number: ES534-BO48440

Expiry date: October 2019

Quality Control Parameters

Lyophilized Powder

MW = 28.5 kDa

Characteristics	Specification	Result
Activity (Dry weight)	≥ 30.0 U/mg	40.8 U/mg
Specific Activity	≥ 40.0 U/mg	48.0 U/mg
RNase	None	
DNase	None	

Authorised by Christopher Weatherall

United Kingdom
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

info.us@bioline.com
Tel: +1 508 880 8990
Fax: +1 508 880 8993

Germany

info.de@bioline.com
Tel: +49 (0)3371 681 229
Fax: +49 (0)3371 681 244

France

info.fr@bioline.com
Tel: +33 (0)1 42 56 04 40
Fax: +33 (0)9 70 06 62 10

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763

Singapore

[Info.sg@bioline.com](mailto:info.sg@bioline.com)
Tel: +65 6774 7196
Fax: +65 6774 6441