



Certificate of Analysis

COA No: CA_BMM-0027

Version: 03

MyTaq Plant-PCR Kit

For Research Use Only

Kit Lot No: BIO-25056_PM843-
BO64120

Storage Conditions: -20°C

Component Lot No: PPX-818410A

Expiry date: November 2020

Quality Control Parameters

Analysis	Specification	Result
Functional	The mix is able to amplify DNA to an acceptable yield from tomato and rice leaves using 40 cycles end-point PCR. Single distinct bands were observed by use of agarose gel electrophoresis (ethidium stained).	Passed
DNase contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation takes 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection 2.5×10^{-3} U DNase.	Passed

Authorised by Christopher Weatherall

United Kingdom
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

info@meridianlifescience.com
Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

info.de@bioline.com
Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

France

info.fr@bioline.com
Tel: +33 (0)1 42 56 04 40
Fax: +33 (0)9 70 06 62 10

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763