

## RNase Inhibitor

Suitable for Research and further Manufacturing Use as an IVD component

Catalog No:	BIO-65028
Lot No:	RA396-B066560
Shipping / Storage Conditions:	-20°C
Component Lot No:	RI-819101A
Expiry date:	February 2021

### Quality Control Parameters

Inhibits a broad spectrum of eukaryotic RNases, including RNases A, B and C to control for contaminants in RT-PCR assays

Analysis	Specification	Result
Inhibition	Test level of inhibition by incubating total RNA with concentration gradient of RNase A. Bands were observed with agarose gel electrophoresis (ethidium stained).	Passed

Authorised by Christopher Weatherall

