

Proteinase K Solution

Suitable for Research and further Manufacturing Use as an IVD component

Catalog No:	BIO-37084
Lot No:	ES537-B068930
Shipping / Storage Conditions:	-20°C
Component Lot No:	PKS-819103A
Expiry date:	April 2021

Quality Control Parameters

RNase and DNase free, ideal for removing endogenous nucleases when purifying native DNA or RNA

Analysis	Specification	Result
Activity	Kinetic colorimetric assay using a substrate to establish kinetics when compared to reference standard.	Passed

Authorised by Christopher Weatherall

