

## JetSeq Flex DNA Library Preparation Kit

Suitable for Research and further Manufacturing Use as an IVD component

Catalog No:	BIO-68027
Lot No:	NG854-B071800
Shipping / Storage Conditions:	-20°C other °
Component Lot No:	NG854-B071800
Expiry date:	July 2021

### Quality Control Parameters

Description	Specification	Result
<p>The JetSeq Flex DNA Library Preparation Kit provides all of the enzymes and reaction buffers for the constructing of next generation sequencing (NGS) libraries.</p> <p>ER Buffer, 5x</p> <p>ER Enzyme Mix</p> <p>Ligation Buffer, 5x</p> <p>Ligase</p> <p>Nuclease Free Water</p> <p>PCR Buffer, 10x</p> <p>DNA Polymerase</p>	<p>JetSeq Flex DNA Library Preparation Kit is used to create a genomic DNA sequencing library. The library was then sequenced on an Illumina sequencer to verify that it meets a series of quality specifications.</p>	Passed
Note to customer	<p>Here provided reagent(s) has(have) been manufactured under the ISO 13485 certified Quality Management System and is(are) suitable for further manufacturing use as an IVD component(s).</p>	N/A

Authorised by Christopher Weatherall

