

**SensiFAST™ Probe Lo-ROX Kit**

Suitable for Research and further Manufacturing Use as an IVD component

Catalog No:	BIO-84005
Lot No:	SF590-B072860
Shipping / Storage Conditions:	-20°C
Component Lot No:	SFPL-819207A
Expiry date:	August 2021

**Quality Control Parameters**

Analysis	Specification	Result
Functional	Quantitative PCR analysis amplifying 6 genes from a dilution series of mouse cDNA under standard conditions. Cq profiles must be consistent for the test and reference sample with ± 0.5 Cq variance.	Passed
DNA contamination	Quantitative PCR analysis with no template. Presence of <i>E. coli</i> and mouse genomic DNA checked. Test sample must amplify in line with control sample.	Passed
Exonuclease contamination	Incubation of a 1Kb double stranded DNA fragment. Incubation for 4 hours at 37°C with dilution series of DNase I. Analysed by agarose gel electrophoresis. Test sample must show less degradation than the limit of detection $2.5 \times 10^{-3}$ U DNase I.	Passed

Authorised by Christopher Weatherall