

ISOLATE II RNA Mini Kit

For Research Use Only

Catalog No:	BIO-52073
Lot No:	IS510-B072900
Shipping / Storage Conditions:	Room Temperature
Component Lot No:	IS510-B072900
Expiry date:	12 months from date of purchase

Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II Filters	n/a	n/a	SPI002517
ISOLATE II RNA Mini Columns	Nucleic acid recovery >90%	Passed	SPB002511
Collection tubes (2 mL)	n/a	n/a	COL001750
Collection tubes (1.5 mL)	n/a	n/a	280119
Lysis Buffer RLY	pH = 6.7-7.3, density	Passed	84826
Wash Buffer RW1	pH = 6.9-7.5, density	Passed	85160
Wash Buffer RW2	pH = 7.5-7.9, density, sterility	Passed	84817
Membrane Desalting Buffer MEM	pH = 7.8-7.1, density	Passed	85172
Reaction Buffer for DNase I RDN	pH = 6.9-7.1, density	Passed	84528
DNase, RNase-free	No RNase activity detectable in connection with our kit	Passed	D1905R
RNase-free Water	Density, no RNase activity, sterility	Passed	84773

Authorised by Christopher Weatherall



United Kingdom
Headquarters UK

info_uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

info@meridianlifescience.com
Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

info.de@bioline.com
Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

France

info.fr@bioline.com
Tel: +33 (0)1 42 56 04 40
Fax: +33 (0)9 70 06 62 10

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763