

Random Hexamer Primer

For Research Use Only

Catalog No:	BIO-38028
Lot No:	RA397-B074920
Shipping / Storage Conditions:	-20°C
Component Lot No:	RHP-919109A
Expiry date:	October 2021

Quality Control Parameters

Analysis	Specification	Result
Functional	A 1Kb fragment is reverse transcribed from Poly A RNA with a dilution series of BioScript™, using standard conditions. Single distinct bands were observed with agarose gel electrophoresis (ethidium stained).	Passed
RNase contamination	Quantitative PCR analysis with high and low RNase standards. Test sample must show less RNase activity than the limit of detection 9.7×10^{-3} ng/μL RNase.	Passed

Authorised by Christopher Weatherall