

100mM dTTP

Shipping: On Dry/Blue Ice

Catalog numbers

Batch No.: See bottle

BIO-39039

1 x 250 µL

Store at -20 °C



Storage and stability:

dTTP is shipped on dry/blue ice. On arrival store at -20 °C for optimum stability. Repeated freeze/thaw cycles should be avoided. For long-term storage, aliquoting is recommended.

Expiry:

When stored under the recommended conditions and handled correctly, quality is retained until the expiry date on the outer box label.

Safety precautions:

Please refer to the material safety data sheet for information.

Quality control specifications:

dTTP is extensively tested for contamination prior to release. dNTPs are >99% pure by HPLC and are free of DNase, RNase, Protease, phosphatase and nicking activity. Lithium salts, >99% deoxynucleoside triphosphates (HPLC, area %), <1% deoxynucleoside monophosphates and deoxynucleoside diphosphates.

Notes:

For research or further manufactured use only.

Trademarks:

HyperLadder is a trademark of Bionline Reagents Ltd.

Description

Ultra-pure dNTPs are enzymatically synthesized from premium quality raw materials, using highly specific production systems in our purpose built facilities. The manufacturing process eliminates impurities and PCR-specific inhibitors such as modified nucleotides, tetraphosphates and pyrophosphates commonly observed in other commercially available dNTP products. dNTPs are purified with quantitative HPLC and possess at least 99% purity.

All dNTPs are supplied as lithium salts in purified water at pH 7.5. Lithium salts have greater resistance to repeated freezing and thawing cycles than sodium salts, and lithium salt dNTP preparations remain sterile over the entire shelf-life due to the bacteriostatic activity of lithium towards various microorganisms.

Features

- Ultra-pure: >99% trisphosphate by HPLC
- Extended shelf-life of 24 months at -20 °C
- Free from PCR inhibitors
- DNase, RNase and Nickase free
- Manufactured in a purpose-built facility
- Custom, bulk and OEM nucleotide service

Applications

- Standard and long range PCR assays
- cDNA synthesis
- qPCR
- Microarrays
- DNA sequencing
- Labelling

dTTP Characteristics

2'-deoxythymidine-5'-triphosphate solution

SPECTRAL & HPLC ANALYSIS	PERFORMANCE
λ_{\max} pH 7.0	267 nm
ϵ at λ_{\max}	9.6 E x mmol ⁻¹ x cm ⁻¹
A_{250}/A_{260}	0.65 ± 0.03
A_{280}/A_{260}	0.73 ± 0.02
Concentration	100 mM ± 5%
Appearance	Clear Colorless Solution
pH of Solution	7.5
dNTP (HPLC Area)	≥ 99%
dNDP (HPLC Area)	< 1%
DNases, RNases, Nicking Activity	Negative

Associated products

Product Name	Pack Size	Catalog No.
dNTP Mix, 100 mM total	1 x 500 µL	BIO-39028
dNTP Set	4 x 250 µL	BIO-39025
HyperLadder™ 1 kb	200 Lanes	BIO-33025
Agarose	500 g	BIO-41025

Product citations

1. Meijer, P.-J., *et al. Meth. Mol. Biol.* **525(3)**, 1-17 (2009).
2. Dellinger, M., *et al. Develop. Biol.* **319(2)**, 309-320 (2008).
3. Hampson, L., *et al. FEBS Lett.* **581(21)**, 3955-3960 (2007).
4. Lloyd, R.E., *et al. Gene.* **172**, 2515-2527 (2006).
5. Tabone, T., *et al. NAR* **34(6)**, e45 (2006).

Bionline Ltd
UNITED KINGDOM

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Meridian Life Science Inc.
USA

Tel: +1 901.382.8716
Fax: +1 901.382.0027

Bionline GmbH
GERMANY

Tel: +49 (0)3371 60222 03
Fax: +49 (0)3371 60222 01

Bionline (Aust) Pty. Ltd
AUSTRALIA

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763