

TriSure™

For research or further manufacturing use only

Catalog No:	BIO-38032
Lot No:	TS-120C
Storage Conditions:	+4°C
Component Lot No:	TS-120C
Expiry date:	12 months from date of purchase

Quality Control Parameters

Analysis	Specification	Result
Functional	RNA extraction from 3T3. rRNA ratio [28s/18s] >1.5 RNA Integrity Number (RIN) >9.0	Passed

QA / QC Representative:



Paul McDermott

Date: 21 April 2021

United Kingdom

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

Australia

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763