

Proteinase K

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-37039
Lot No:	ES534-B096820
Storage Conditions:	-20°C
Component Lot No:	ES534-B096820
Expiry date:	August 2023

Quality Control Parameters

Lyophilized Powder

MW = 28.5 kDa

Characteristics	Specification	Result
Activity (Dry weight)	≥ 30.0 U/mg	30.56 U/mg
Specific Activity	≥ 40.0 U/mgP	43.30 U/mgP
RNase	None	
DNase	None	

QA / QC Representative:



Andrew Galeeba-M

 Date: 12th August 2021

United Kingdom

 Tel: +44 (0)20 8830 5300
 Fax: +44 (0)20 8452 2822

USA

 Tel: +1 901.382.8716
 Fax: +1 901.382.0027

Germany

 Tel: +49 (0)3371 60222 00
 Fax: +49 (0)3371 60222 01

Australia

 Tel: +61 (0)2 9209 4180
 Fax: +61 (0)2 9209 4763