

**Proteinase K**

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-37039
Lot No:	ES534-B125370
Storage Conditions:	-20°C
Component Lot No:	ES534-B125370
Expiry date:	April 2026


**Quality Control Parameters**

**Lyophilized Powder**

**MW = 28.5 kDa**

Characteristics	Specification	Result
Activity (Dry weight)	≥ 30.0 U/mg	45.2 U/mg
Specific Activity	≥ 40.0 U/mgP	≥ 40.0 U/mgP
RNase	None	Passed
DNase	None	Passed

QA / QC Representative:



J. Rahnenführer

Date: 9<sup>th</sup> February 2024

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01