

**RNase Inhibitor**

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-65028
Lot No:	RA396-B105200
Storage Conditions:	-20°C
Component Lot No:	RI-122302A
Expiry date:	March 2024

**Quality Control Parameters**

Inhibits a broad spectrum of eukaryotic RNases, including RNases A, B and C to control for contaminants in RT-PCR assays

Analysis	Specification	Result
Inhibition	Test level of inhibition by incubating total RNA with concentration gradient of RNase A. Bands were observed with agarose gel electrophoresis (ethidium stained).	Passed

QA / QC Representative:



Alberta Newton

Date: 24<sup>th</sup> February 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763