

**RNase Inhibitor**

Suitable for Research and further Manufacturing Use


Catalog No:	BIO-65028
Lot No:	RA396-B125760
Storage Conditions:	-20°C
Component Lot No:	RI-124302A
Expiry date:	March 2026

**Quality Control Parameters**

Inhibits a broad spectrum of eukaryotic RNases, including RNases A, B and C to control for contaminants in RT-PCR assays

Analysis	Specification	Result
Inhibition	Test level of inhibition by incubating total RNA with concentration gradient of RNase A. Bands were observed with agarose gel electrophoresis (ethidium stained).	Passed

QA / QC Representative:



J. Rahnenführer

Date: 29<sup>th</sup> February 2024

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01