

# ISOLATE II Gel and PCR Kit


Suitable for Research and further Manufacturing Use

Catalog No:	BIO-52059
Lot No:	B123760
Storage Conditions:	Room Temperature
Component Lot No:	IS496-B123760
Expiry date:	12 month from date of purchase

## Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II PCR and Gel Columns	Nucleic acid recovery >90%	Passed	SPI003275
Collection tubes	n/a	n/a	COL002027
Binding Buffer CB	pH = 4.9-5.1, density	Passed	115654
Wash Buffer CW (conc)	pH = 7.5-7.9, density	Passed	115313
Elution Buffer C	pH = 8.3-8.7, density	Passed	115529

QA / QC Representative:



J. Rahnenführer

Date: 21<sup>st</sup> March 2024

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01