

## ISOLATE II Gel and PCR Kit

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-52059
Lot No:	IS496-B100130
Storage Conditions:	Room Temperature
Component Lot No:	IS496-B100130
Expiry date:	12 months from date of purchase

### Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II PCR and Gel Columns	Nucleic acid recovery >90%	Passed	SPI002829
Collection tubes	n/a	n/a	COL001895
Binding Buffer CB	pH = 4.8-5.2, density	Passed	102157
Wash Buffer CW (conc)	pH = 7.5-7.9, density	Passed	100892
Elution Buffer C	pH = 8.3-8.7, density	Passed	101807

QA / QC Representative:



Andrew Galeeba-M

Date: 6<sup>th</sup> December 2021

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763