

ISOLATE II Gel and PCR Kit

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-52060
Lot No:	IS497-B120510
Storage Conditions:	Room Temperature
Component Lot No:	IS497-B120510
Expiry date:	12 months from date of purchase

Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II PCR and Gel Columns	Nucleic acid recovery >90%	Passed	SPI003206
Collection tubes	n/a	n/a	COL001967
Binding Buffer CB	pH = 4.9-5.1, density	Passed	112870
Wash Buffer CW (conc)	pH = 7.5-7.9, density	Passed	112697
Elution Buffer C	pH = 8.3-8.7, density	Passed	112937

QA / QC Representative:



Andrew Galeeba-M

 Date: 15th November 2023

United Kingdom

 Tel: +44 (0)20 8830 5300
 Fax: +44 (0)20 8452 2822

USA

 Tel: +1 901.382.8716
 Fax: +1 901.382.0027

Germany

 Tel: +49 (0)3371 60222 00
 Fax: +49 (0)3371 60222 01