

## ISOLATE II RNA Mini Kit

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-52072
Lot No:	IS509-B091140
Storage Conditions:	Room Temperature
Component Lot No:	IS509-B091140
Expiry date:	12 months from date of purchase

### Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II Filters	n/a	n/a	SPI002748
ISOLATE II RNA Mini Columns	Nucleic acid recovery >90%	Passed	SPB002716
Collection tubes (2 mL)	n/a	n/a	COL001864
Collection tubes (1.5 mL)	n/a	n/a	191020
Lysis Buffer RLY	pH = 6.7-7.3, density	Passed	95221
Wash Buffer RW1	pH = 6.9-7.5, density	Passed	95213
Wash Buffer RW2	pH = 7.5-7.9, density, sterility	Passed	95543
Membrane Desalting Buffer MEM	pH = 7.8-7.1, density	Passed	95189
Reaction Buffer for DNase I RDN	pH = 6.9-7.1, density	Passed	94786
DNase, RNase-free	No RNase activity detectable in connection with our kit	Passed	D2013R
RNase-free Water	Density, no RNase activity, sterility	Passed	95272

QA / QC Representative:



Andrew Galeeba-M

Date: 17<sup>th</sup> February 2021

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763