

JetSeq Clean

Suitable for Research and further Manufacturing Use

Catalog No:	BIO-68032
Lot No:	NG876-B096800
Storage Conditions:	at +4°C
Component Lot No:	NG876-B096800
Expiry date:	7 th June 2023

Quality Control Parameters

Analysis	Specification	Result
Double-sided size selection	JetSeq Clean beads are able to select a specific size range from a mixture of DNA sizes, in standard double-sided size selection conditions.	PASS

QA / QC Representative:



Andrew Galeeba-M

 Date: 28th June 2021

United Kingdom

 Tel: +44 (0)20 8830 5300
 Fax: +44 (0)20 8452 2822

USA

 Tel: +1 901.382.8716
 Fax: +1 901.382.0027

Germany

 Tel: +49 (0)3371 60222 00
 Fax: +49 (0)3371 60222 01

Australia

 Tel: +61 (0)2 9209 4180
 Fax: +61 (0)2 9209 4763