

JetSeq™ Clean

Suitable for Research and further Manufacturing Use

| | |
|---------------------|-------------------------------|
| Catalog No: | BIO-68032 |
| Lot No: | NG876-B104860 |
| Storage Conditions: | at +4°C |
| Component Lot No: | NG876-B104860 |
| Expiry date: | 19 th January 2024 |

Quality Control Parameters

| Analysis | Specification | Result |
|-----------------------------|--|-------------|
| Double-sided size selection | JetSeq Clean beads are able to select a specific size range from a mixture of DNA sizes, in standard double-sided size selection conditions. | PASS |

QA / QC Representative:



Andrew Galeba-M

Date: 8th February 2022

United Kingdom

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

Australia

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763