

**IPTG**

Suitable for Research and further Manufacturing Use

Catalog No.:	BIO-37036
Lot No	ES532-B108210
Storage Conditions:	-20°C
Component Lot No:	ES532-B108210
Expiry date:	February 2024

**Quality Control Parameters**

Analysis	Specification	Result
Appearance	White powder	Passed
Purity	>99.5%	Passed
Melting range	111°C ± 2°C	Passed
Optical rotation	-30-32	Passed

QA / QC Representative:



Andrew Galeeba-M

Date: 9<sup>th</sup> June 2022

**United Kingdom**

Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
Fax: +49 (0)3371 60222 01

**Australia**

Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763