

## Proteinase K Solution

For research or further manufacturing use only


Catalog No:	BIO-37085
Lot No:	ES538-B126440
Storage Conditions:	-20°C
Component Lot No:	PKS-324104A
Expiry date:	May 2026

### Quality Control Parameters

**Composition:** 20 mg/ml Proteinase K, 10 mM Tris/Cl pH 7.5, 2 mM CaCl<sub>2</sub>, 50% Glycerol

Analysis	Specification	Result
Activity	Kinetic colorimetric assay using a substrate to establish kinetics when compared to reference standard.	Passed

QA / QC Representative:



J. Rahnenführer

 Date: 9<sup>th</sup> April 2024

**United Kingdom**

Tel: +44 (0)20 8830 5300  
 Fax: +44 (0)20 8452 2822

**USA**

Tel: +1 901.382.8716  
 Fax: +1 901.382.0027

**Germany**

Tel: +49 (0)3371 60222 00  
 Fax: +49 (0)3371 60222 01