

Proteinase K Solution

For research or further manufacturing use only


Catalog No:	BIO-37084
Lot No:	ES537-B106850
Storage Conditions:	-20°C
Component Lot No:	PKS-022104A
Expiry date:	May 2024

Composition: 20 mg/ml Proteinase K, 10 mM Tris/Cl pH 7.5, 2 mM CaCl₂, 50% Glycerol

Quality Control Parameters

Analysis	Specification	Result
Activity	Kinetic colorimetric assay using a substrate to establish kinetics when compared to reference standard.	Passed

QA / QC Representative:



Alberta Newton

Date: 14th April 2022

United Kingdom

Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

USA

Tel: +1 901.382.8716
Fax: +1 901.382.0027

Germany

Tel: +49 (0)3371 60222 00
Fax: +49 (0)3371 60222 01

Australia

Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763