

## Proteinase K Solution

For research or further manufacturing use only

Catalog No:	BIO-37084
Lot No:	ES537-B115090
Storage Conditions:	-20°C
Component Lot No:	PKS-123202A
Expiry date:	March 2025

**Composition:** 20 mg/ml Proteinase K, 10 mM Tris/Cl pH 7.5, 2 mM CaCl<sub>2</sub>, 50% Glycerol

### Quality Control Parameters

Analysis	Specification	Result
Activity	Kinetic colorimetric assay using a substrate to establish kinetics when compared to reference standard.	Passed

QA / QC Representative:



Andrew Galeeba-M

 Date: 28<sup>th</sup> February 2023

United Kingdom

 Tel: +44 (0)20 8830 5300  
 Fax: +44 (0)20 8452 2822

USA

 Tel: +1 901.382.8716  
 Fax: +1 901.382.0027

Germany

 Tel: +49 (0)3371 60222 00  
 Fax: +49 (0)3371 60222 01

Australia

 Tel: +61 (0)2 9209 4180  
 Fax: +61 (0)2 9209 4763