

Genotyping Blood qPCR Mix

Cat No: ODX001

Shipping conditions: -20 °C

Lot number: RD123F

Test Parameters

Analysis	Specification	Result
Functionality	<p>Quantitative PCR analysis amplifying target from a dilution series of template under standard conditions.</p> <p>Relative activity is measured against a reference by a functional qPCR assay.</p> <p>Test Criteria: Test and reference amplification profiles must align; ΔCt (Test-Reference) < 1</p>	Passed

Authorised by Michele Amasio
Sr. Manager, R&D

Michele Amasio 28/05/2023