



Certificate of Analysis

COA No: CA_RES-0021

Version: 01

ISOLATE II Gel and PCR Kit

For Research Use Only

Storage
Conditions:

Room Temperature

Lot number:

IS496-BO25550

Expiry date:

12 months from date of
purchase

Quality Control Parameters

| Component | Specification | Result | Batch |
|--------------------------------|----------------------------|--------|----------|
| ISOLATE II PCR and Gel Columns | Nucleic acid recovery >90% | Passed | SC-115D |
| Collection tubes | n/a | n/a | OT-115D |
| Binding Buffer CB | pH = 4.8-5.2, density | Passed | CB-115D |
| Wash Buffer CW (conc) | pH = 7.5-7.9, density | Passed | CW-115D |
| Elution Buffer C | pH = 8.3-8.7, density | Passed | EBC-115D |

Authorised by Jade James

Europe
Headquarters UK

info.uk@bioline.com
Tel: +44 (0)20 8830 5300
Fax: +44 (0)20 8452 2822

Europe
Germany

info.de@bioline.com
Tel: +49 (0)3371 681 229
Fax: +49 (0)3371 681 244

America

info.us@bioline.com
Tel: +1 508 880 8990
Fax: +1 508 880 8993

Australia

info.aust@bioline.com
Tel: +61 (0)2 9209 4180
Fax: +61 (0)2 9209 4763