



## Certificate of Analysis

COA No: CA\_RES-0020

Version: 01

### ISOLATE II Gel and PCR Kit

For Research Use Only

Storage  
Conditions:

Room Temperature

Lot number:

IS495-BO26400

Expiry date:

12 months from date of  
purchase

### Quality Control Parameters

Component	Specification	Result	Batch
ISOLATE II PCR and Gel Columns	Nucleic acid recovery >90%	Passed	SC-115F
Collection tubes	n/a	n/a	OT-115F
Binding Buffer CB	pH = 4.8-5.2, density	Passed	CB-115F
Wash Buffer CW (conc)	pH = 7.5-7.9, density	Passed	CW-115F
Elution Buffer C	pH = 8.3-8.7, density	Passed	EBC-115F

Authorised by Jade James

**Europe**  
Headquarters UK

[info.uk@bioline.com](mailto:info.uk@bioline.com)  
Tel: +44 (0)20 8830 5300  
Fax: +44 (0)20 8452 2822

**Europe**  
Germany

[info.de@bioline.com](mailto:info.de@bioline.com)  
Tel: +49 (0)3371 681 229  
Fax: +49 (0)3371 681 244

**America**

[info.us@bioline.com](mailto:info.us@bioline.com)  
Tel: +1 508 880 8990  
Fax: +1 508 880 8993

**Australia**

[info.aust@bioline.com](mailto:info.aust@bioline.com)  
Tel: +61 (0)2 9209 4180  
Fax: +61 (0)2 9209 4763